

# CAREPLUS

YOUR ONE-STOP CENTER FOR QUALITY MALAYSIAN MADE PROTECTION PRODUCTS

## Product Information Service

### NITRILE EXAMINATION GLOVES—POWDERED

(Product Code: NE01)

CAREPLUS (M) SDN BHD (Company No. 212677-K)

Lot 110, Lorong Senawang 4/3, Off Jalan Senawang Empat, Senawang Industrial Estate, 70450 Seremban, N. Sembilan, Malaysia.

Tel: (60) 6 6772781 Fax: (60) 6 6772780 Email: careplus@po.jaring.my Website: www.careplus.com.my

**Reference International Standards:** ASTM D3578-00a, EN455-1:2000, EN455-2:2000 and EN455-3:1999.

**Sampling Plan:** MIL STD 105E or ISO 2859

**Quality System:** ISO9002:1994 and EN46002:1996

**Other Quality Endorsements:** CE Mark

**Product Description:** Ambidextrous, Single-Use, Powdered, Beaded Cuff, Micro or Finger-Textured, White or Blue Color.

**Base Material:** 100% Acrylonitrile-Butadiene Synthetic Latex.

**Glove Length:** 240mm minimum  
(also available in 300mm)

**Typical Thickness (Single wall):** 0.11mm at cuff, 0.14mm at palm and 0.16mm at finger  $\pm$  0.01 mm.

	<u>AQL Limit</u>	<u>Inspection Level</u>
Pinholes	1.5	G1
Visual Defects—Major	2.5	G1
Visual Defects—Minor	4.0	G1
Physical Dimensions	4.0	S2
Physical Properties	4.0	S2
Quantity Count	4.0	S1

(Note: gloves are packed to 100 pieces per box as a standard practice)

**Powder Content:** Maximum 120 mg per glove of USP absorbable dusting powder

**Typical Physical Properties**

	<u>Before Ageing</u>	<u>After Ageing</u>
--	----------------------	---------------------

Tensile Strength	14 Mpa	14 Mpa
------------------	--------	--------

Ultimate Elongation	700 %	500 %
---------------------	-------	-------

		<u>Width</u>
Typical Glove Width	X-Small	70—80mm
	Small	80—90mm
	Medium	90—100mm
	Large	105—115mm
	X-Large	110—120mm

**Packing:** 100 gloves/box, 10 boxes/ carton, 1200 /2480 cartons per 20'/40' container.

